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**Indications:**

For the therapy of aches and pains with heat or cold. Fast and easy application.

**Important:**

- 1- Follow instructions for application.
- 2- Do not use gel packs if punctured or torn.
- 3- Persons with circulatory impairment please consult their physician before using hot/cold packs.
- 4- FLEXXUM® gel packs are non-toxic, but should not be eaten. By ingestion, rinse mouth cavity. Induce vomiting when larger quantities ingested. Supervise use with small children.

**Cold use:**

Store FLEXXUM® gel packs flat in your freezer or refrigerator. The gel packs will stay pliable to -25°C (for FLEXXUM®) or to -12°C (for FLEXXUM® B).

**Hot use:**

**Heating in hot water:** Boil enough hot water to cover FLEXXUM® gel packs fully. Then take the boiling water off the heat and place the gel pack carefully into the hot water. Heat approx. 7 min. (10 min. if previously frozen).

**Heating in microwave:** FLEXXUM® may be heated in microwave. Heating time depends on the make of your oven and the size of the gel pack. Further instructions on the back.

Most microwave ovens come in 350 to 900 watt strength, but some go higher. Because of this large variation, heating times can differ. Thus, it is very important to be present during the entire first heating process with your microwave oven. This will allow you to prolong or shorten the heating time as needed.

With the first heating by microwave, always start with the shorter of the heating times given below. If this leads to insufficient heating, heat in additional steps of 10-30 sec. (according to size of the gel pack) until the maximum time given. Then add up the total heating time in your microwave and use this time from now on. This time should stay the same for future uses on the same microwave.

Art. No.	Size	Heating times (in seconds)
040	13x15cm	20 - 40
041	13x30cm	60 - 120
043	19x30cm	120 - 240
044	30x38cm	180 - 360
045/057	7x38cm	30 - 60
047/053	7x19cm	15 - 30
048	26x30cm	120 - 240
049	7x9,5cm	10 - 20
050	42x38cm	120 - 240
051	84x38cm	no microwave
054	19x42cm	90 - 180
052	15x19cm	45 - 90
058	19x45cm	150 - 300
060	15x38cm	90 - 180
063	7x25cm	20 - 40
064	7x13cm	10 - 20
066	10x26cm	40 - 80

**Attention:** If the plastic foil of the gel starts to rise, the content will be too hot, i.e. near boiling. Turn off microwave right away and handle gel pack very carefully. When warming next time, use 10-30 sec. shorter heating time.

**Please consider** that products heated by microwave will increase in temperature within the following few minutes. It is therefore best to wait a few minutes before heating some more.

**Re-heating of a tepid pack:** Half of regular heating time.

**>>>> Caution! Very hot or very cold gel packs may lead to burns. Never place gel packs directly onto bare skin.**

OVERHEATING WILL INVALIDATE OUR PRODUCT GUARANTEE

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